



Reference: PTXBC063489

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	FP2653; transmembrane and ubiquitin-like domain-containing protein 2; transmembrane and ubiquitin like domain containing 2		

For research use only.

GENE DESCRIPTION	transmembrane and ubiquitin-like domain containing 2
GENE ID	79089
ORIGIN SPECIES	Human
ACCESSIONS	BC063489
VECTOR	pDONR223
NCBI SYMBOL	TMUB2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]