

FAM122C-family with sequence similarity 122C Gene

Reference: PTXBC063491

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein FAM122C; family with sequence similarity 122C		
For research use only.			

GENE DESCRIPTION	family with sequence similarity 122C
GENE ID	159091
ORIGIN SPECIES	Human
ACCESSIONS	BC063491
VECTOR	pDONR223
NCBI SYMBOL	FAM122C
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcacaggagaaaatgaaactaggttcaagtcgctgccgagttccactaccgcagacggcaacattctgagaagagtcacagtgccctttgatcaatggacttggttttaattcacaggtgttgcaagctgacatgtaagaattaggacaaacagaacaacatttaggaatcgacgtctctgtttgcccactctccctttcatggttccatcagccgcttcatcaaatcaaacaggagaagccatggatttaataaataagagaacaatgtctgaatggaagctacaaagttagatacagataagtcactcttgggaagaaggcttgaactggtgaaatggcatttcaatatcaaccaaaaagattttccaaggcacaaccaacatgctttctctgatacttcccaactgtcagaaaaatgtgtactactccagctacttttgatggaacgacagcaatgctggatcttctggaattcgtcagctgaaatcggcactgtcacaaactctctgtgtcaccttctgatactggttctcattgttctagtagacaaacttcaactaa