

BCOR-BCL6 co-repressor Gene

Reference: PTXBC063536

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ANOP2; MAA2; MCOPS2; BCL-6 corepressor; BCL-6 interacting corepressor; BCL6 corepressor		

For research use only.

GENE DESCRIPTION	BCL6 co-repressor
GENE ID	54880
ORIGIN SPECIES	Human
ACCESSIONS	BC063536
VECTOR	pDONR223
NCBI SYMBOL	BCOR
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacccacagtgaacttatgaaaagttcttaacagattattaaatgacctccagggtcgcaatgatgatgacgccagtgccacttgggacttctatggcagctctgtttgtgaaccagatgatgaaagtggctatgatgttttagccaacccccaggaccagaagaccaggatgatga
cgatgcctatagcagatgtgttgaatttgaatttcagagacccccctctaccgtgttataacatccaagtatctgtgctcaggggccacgaaactggctactgtttcggatgctcttaagaaatgaaaatgtcccccataatttcgtgcaattttcacaacgtggaattgtcaccattgca
gaggcagaattttatcggcaggtttctgcaagtccttctgtctcctgctccaaagacctggaagccttcaaccctgaaagtaaggagctgttagatctgttggaattcacgaacgaaatcagactctgtgggtcctctgtagagtggtcccccagtgatctggcctcagacaactactgg
tga