



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

PAX3-paired box 3 Gene

Reference: PTXBC063547

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CDHS; HUP2; WS1; WS3; paired box protein Pax-3; paired box homeotic gene 3; paired domain gene 3; paired domain gene HuP2; paired box 3		

For research use only.

GENE DESCRIPTION	paired box 3
GENE ID	5077
ORIGIN SPECIES	Human
ACCESSIONS	BC063547
VECTOR	pDONR223
NCBI SYMBOL	PAX3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]