

UPK1B-uroplakin 1B Gene

Reference: PTXBC063568

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	TSPAN20; UPIB; UPK1; uroplakin-1b; UP1b; tetraspan; tetraspanin-20; tspan-20; uroplakin Ib; uroplakin 1B		

For research use only.

GENE DESCRIPTION	uroplakin 1B
GENE ID	7348
ORIGIN SPECIES	Human
ACCESSIONS	BC063568
VECTOR	pDONR223
NCBI SYMBOL	UPK1B
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccaaagacaactcaactgttcgttcttcaggcctgctgatttttgaaatgtgattattggtgttcggcattgccctgactcggagtgcatctctttgtatctgaccaacacagcctctaccactgcttgaagccaccgacaacgatgacatctatgggctgctggatcgcca
tatttggggcatctgctcttctgctgtctgttctagcattgtaggcacatgaagtcagcaggaataattctctgctgatttcatctgatgtttatagtatatgctttgaagtggtcatcttgcacagcagcaacacacagacattttcacaccaacctcttctgaagcagatgctag
agaggtaccaaacaacagccctccaaacaatgatgaccagtggaaaaacaatggagtcacaaaacctgggacaggctcatgctccaggacaattgctgtggcgtaaatgggtccatcagactggcaaaatacacatctgcttccggactgagaataatgatgctgactatccctggcc
tcgtcaatgctgtgtatgaacaatcttaagaacctctcaacctggaggctgtgaaactaggcgtgctgttttatcacaatcagggtgctatgaactgatctctggccaatgaaccgacacgcctggggggttgcctgtgttgattgccaattctctgctggaattttgggtctcctgggt
accatgttctactggagcagaattgaattataa