

PASK-PAS domain containing serine/threonine kinase Gene

Reference: PTXBC063585

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PASKIN; STK37; PAS domain-containing serine/threonine-protein kinase; per-arnt-sim (PAS) domain kinase; PAS domain containing serine/threonine kinase		

For research use only.

GENE DESCRIPTION	PAS domain containing serine/threonine kinase
GENE ID	23178
ORIGIN SPECIES	Human
ACCESSIONS	BC063585
VECTOR	pDONR223
NCBI SYMBOL	PASK
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1