

PTPRK-protein tyrosine phosphatase, receptor type, K Gene

Reference: PTXBC063596

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	R-PTP-kappa; receptor-type tyrosine-protein phosphatase kappa; dJ480J14.2.1 (protein tyrosine phosphatase, receptor type, K (R-PTP-KAPPA, protein tyrosine phosphatase kappa , protein tyrosine phosphatase kappa; protein-tyrosine phosphatase kappa; protein-		

For research use only.

GENE DESCRIPTION	protein tyrosine phosphatase, receptor type, K
GENE ID	5796
ORIGIN SPECIES	Human
ACCESSIONS	BC063596
VECTOR	pDONR223
NCBI SYMBOL	PTPRK
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggatacgaactgcggcggcgctgctgtttgtggcgtcttgctctctctcttggcctctctgggatcgcccaaggccaggtctccgcagtgctgtgctgtgtgctgtaactttgctgctgtgtgtctgctgtctgccattccattcgccatctatcacgcacagaaacct
ggccagatcatgaaggcgattgaaaagggtgctgtactttga