

HSPBAP1-HSPB (heat shock 27kDa) associated protein 1 Gene

Reference: PTXBC063629

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PASS1; HSPB1-associated protein 1; 27 kDa heat shock protein-associated protein 1; HSPB (heat shock 27kDa) associated protein 1; protein associating with small stress protein PASS1; HSPB1 associated protein 1		

For research use only.

GENE DESCRIPTION	HSPB (heat shock 27kDa) associated protein 1
GENE ID	79663
ORIGIN SPECIES	Human
ACCESSIONS	BC063629
VECTOR	pDONR223
NCBI SYMBOL	HSPBAP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagcaggtctccgagcgaccactcctgtgatcgttgcggctggggctggaggagggaaggtgaacatgtcaaaccttttaagccagagaagcaaaagaaattatcatgtctttacaacaacctgcaatctctgaacatgggtgttgattggccagcagcagactggaatgctaa
atacctttcgcaggtccttcctcatggcaagcagatacgaatgggafgaaaagcatgagcacagttcctcagtttgaactacatgtaattacgtagaagctacactcgaagatttctgacctggaactgtgaccagtctagtatttctggaccatttagagattatgaccattccaagttct
ggccttatgtgactataaattttgtcagttatftgaagacaagacagatctttccaggatgtgaatggctgacttcgggttctctggaagaaatggacaggaagatgattgtgattggctcttgggagccacacacctgtcatctggactcctatggttgaactgggtattccagg
tacaaggaagcctggagtgcaatggcatgatcatagctcctggacctcaagccatcctcctgcgcagcctctcaagtatctaggactacaggaacgatgcatctcttctcctga