



Reference: PTXBC063633

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PDXDC2; PDXDC2P; pyridoxal-dependent decarboxylase domain containing 2; pyridoxal-dependent decarboxylase domain-containing protein 1		

For research use only.

GENE DESCRIPTION	pyridoxal-dependent decarboxylase domain containing 2
GENE ID	283970
ORIGIN SPECIES	Human
ACCESSIONS	BC063633
VECTOR	pDONR223
NCBI SYMBOL	PDXDC2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]