

CCM2-cerebral cavernous malformation 2 Gene

Reference: PTXBC063663

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CCM2 scaffolding protein; C7orf22; OSM; PP10187; cerebral cavernous malformations 2 protein; cerebral cavernous malformation 2; malcavernin; osmosensing scaffold for MEKK3		

For research use only.

GENE DESCRIPTION	cerebral cavernous malformation 2
GENE ID	83605
ORIGIN SPECIES	Human
ACCESSIONS	BC063663
VECTOR	pDONR223
NCBI SYMBOL	CCM2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaagaggaggcaagaaggcaagaagcctggaattgtctgcatttaaacgagtagttctctaaagggtgaaaagagtagagataagaagcccatgagaaggtgacagagaggcgcctctgcacactgtggtgtgtcattgcctgagcgcgtcgagccagacagactgctgagcactatattgagaaggaggtaaaggctccattaa