

CREB3L2-cAMP responsive element binding protein 3-like 2 Gene

Reference: PTXBC063666

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	cyclic AMP-responsive element-binding protein 3-like protein 2; B-ZIB transcription factor; BBF2 human homolog on chromosome 7; FUS/BBF2H7 protein; TCAG_1951439; basic transcription factor 2; cAMP-responsive element-binding protein 3-like protein 2; cAMP		

For research use only.

GENE DESCRIPTION	cAMP responsive element binding protein 3-like 2
GENE ID	64764
ORIGIN SPECIES	Human
ACCESSIONS	BC063666
VECTOR	pDONR223
NCBI SYMBOL	CREB3L2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggtgctggagagcggggagcagggcgctgctgcagtgaggaccgcaagctgagcgagctgtcagagccggggacggcgaggccctcatgtaccacacgcactctcagaactctggatgagtttccagaacgtctgggtcagctcctgaatgatcctttctcagagaa
gagtgctcaatggagtggaaccttccccgacgtccccggcgctctcatccaggctgagcagactactcctgtgagggagcctcgggcccagtcgccccttaccacattaccagtgacagcttcaatgacgatgagtggaagtggagaaatggtacgtgtacagactcct
tcaacatcatcaagacagagccaggttacagacgaaccaccccgactgttccgtgtgctcctgacctgacctacagccatctccacccgttgaaaaaggaggaaacctcctctggaaatgaacactggggtgattcctgtgcccagaccattatcctaaaaatgaagctggagcctcat
gaagtggatcagtttctaaacttctcctaaagaaggtgtgtgtgccccctgtgtcccttgggttatgatattggtctctgggtctacagagagggaatatggcgagagagctgggatgagttgtaccacagatgtttagctggctttatgaaatagctctgttctaaaaataaaatttg
cttccaaataa