



Reference: PTXBC063676

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ELK4, ETS transcription factor; ELK4, ETS-domain protein (SRF accessory protein 1); SAP1; ETS domain-containing protein Elk-4; SAP-1; serum response factor accessory protein 1		

For research use only.

GENE DESCRIPTION	ELK4, ETS-domain protein (SRF accessory protein 1)
GENE ID	2005
ORIGIN SPECIES	Human
ACCESSIONS	BC063676
VECTOR	pDONR223
NCBI SYMBOL	ELK4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]