

LOC595101-PI-3-kinase-related kinase SMG-1 pseudogene Gene

Reference: PTXBC063703

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms			
For research use only.			

GENE DESCRIPTION	PI-3-kinase-related kinase SMG-1 pseudogene
GENE ID	595101
ORIGIN SPECIES	Human
ACCESSIONS	BC063703
VECTOR	pDONR223
NCBI SYMBOL	LOC595101
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgagccgcagagccccgggtctggctgagcagcggcgccacacactatcgggagctggaatgactggcaaccagaaactgatagtcacacagctgacccaggttaattaaaatattcttcacagagatagaggtggtttctctctacggactgcaaccttcaaatcagctgtg
gtgtctcggaagggcagcagatgataccagagtcacgctgacatacagaatgacgaaaagggtggctacagtgtaatggagatctggggaataacttatgtcgggaagtcgttggggcaagagctgagggttaacaatgtgaccagccctgagttcaccaggtgtcagcatggcagt
cgtgcttagccaccaagacatgaggaatcacaggagagatcgtcttattgtgatgagtcctgactgtcaaatctcttcggagatcacccgggaagacgacagagactga