

TMEM91-transmembrane protein 91 Gene

Reference: PTXBC063705

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DSPC3; IFITMD6; transmembrane protein 91; dispanin subfamily C member 3; interferon induced transmembrane protein domain containing 6		

For research use only.

GENE DESCRIPTION	transmembrane protein 91
GENE ID	641649
ORIGIN SPECIES	Human
ACCESSIONS	BC063705
VECTOR	pDONR223
NCBI SYMBOL	TMEM91
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacagccctagcttctgtagcttcaacagcctctgctggaggccacagaatgtgagaccctgccagaagcctggcagcagctgagctgggtccccccttaagagagatagccttgcgagtcctgaggggtttgagttcctgtcaccgcctctccctccgtgagcgtggcct
gggggaaccaaggcccccctgagttgaggacatgtcatccagtgacagtgactcggactgggatggaggcagccgtcttccacattctaccccacgaccacctcggcttgctgtcttccatgctgtgtgtttctgcccgttggcatcgtgcttctgtctagcccagaagcatcccc
agtgcctgacctatagtaggtgctacgtatctgcttttgggaatacatgcatga