



Reference: PTXBC063795

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ERR1; ERRa; ERRalpha; ESRL1; NR3B1; steroid hormone receptor ERR1; ERR-alpha; estrogen receptor-like 1; estrogen-related nuclear receptor alpha; nuclear receptor subfamily 3 group B member 1; estrogen related receptor alpha		

For research use only.

GENE DESCRIPTION	estrogen-related receptor alpha
GENE ID	2101
ORIGIN SPECIES	Human
ACCESSIONS	BC063795
VECTOR	pDONR223
NCBI SYMBOL	ESRRA
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]