

DNAJC24-DnaJ (Hsp40) homolog, subfamily C, member 24 Gene

Reference: PTXBC063804

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	DPH4; JJJ3; ZCSL3; dnaJ homolog subfamily C member 24; 1700030A21Rik; CSL-type zinc finger-containing protein 3; DPH4 homolog (JJJ3, S. cerevisiae); DPH4, JJJ3 homolog; DnaJ (Hsp40) homolog, subfamily C, member 24; diphthamide biosynthesis protein 4; zinc		

For research use only.

GENE DESCRIPTION	DnaJ (Hsp40) homolog, subfamily C, member 24
GENE ID	120526
ORIGIN SPECIES	Human
ACCESSIONS	BC063804
VECTOR	pDONR223
NCBI SYMBOL	DNAJC24
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcgggtgagcagatgccaaaaaggattggtacagcatcctgggagcagaccatctgcaatatatcagacctaaacaaaaatcaaaaactatattaatgtatcatccagataaacaagtacagatgtaccagcaggaacagtggaggaatgtgtacagaagttcatcgaaatt
gatcaagcatggaaaattctaggaaatgaagagacaaaaagagatgatgacctgcagcgggtgtgaagatgatctaagaatgtaggaccagtagatgctcaagtatatcttgaagaatgtcttggatgaaggtgatcactcttttatctgagttgcagatgtgtggaaaatcacagtgtttc
caaggatgaagcggagaagtttagcctgatttctgtgatacatgttcactaattatagaactcttcattataactaa