



Reference: PTXBC063815

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	INT3; neurogenic locus notch homolog protein 4; Notch homolog 4; notch 4		
For research use only.			

GENE DESCRIPTION	Notch homolog 4 (Drosophila)
GENE ID	4855
ORIGIN SPECIES	Human
ACCESSIONS	BC063815
VECTOR	pDONR223
NCBI SYMBOL	NOTCH4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]