



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

GBP3-guanylate binding protein 3 Gene

Reference: PTXBC063819

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	guanylate-binding protein 3; GTP-binding protein 3; guanine nucleotide-binding protein 3; guanylate-binding protein 3 delta C; guanylate binding protein 3		

For research use only.

GENE DESCRIPTION	guanylate binding protein 3
GENE ID	2635
ORIGIN SPECIES	Human
ACCESSIONS	BC063819
VECTOR	pDONR223
NCBI SYMBOL	GBP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgtgcctcagagatcacatgtacagcgcccaattgtgctctatgagacaactaaagggaactgtgtggcaatccagaagctctgaaataacctgtctgcattacacagcctgtgtgtgtgtggcaattgtggcgctctaccgcacaggaataactactgatgaacaagctagctgtgga
agaataaaggctctctctgtggctccacagtgaactctcacaccaaaaggaaactgtgatgtgtgtgtgctctccccaaaaagccagaacacaccttagtctctgtgactctgaggcgctgggagatgaagaaggggtgacaaccagaatgactctggatctcacctgtggcgtctctc
ctgacagcactctctctgtgataaagtatgaaacataacacagcaggtatggacaacactgtactatctgacagagctgacacactgaaatccgatcaaaatctctacctgtatgaaagatgaaagttaggaactgacttgtgagctctctccagatattgtgtggacactgagagattt
ctctcctggacttggagcagatggacaaccctccacacagatgacttacctggaattatctccgaagtcaacgcaggtaccacgtcaaaaataaaaataaaatttaactctgcccgcactctgtatccggaaacttctcccaagaataaaatgttttgtctgactgtcccatccaccgcaggaagc

More info: Please contact us for any question: contact@proteogenix.fr