



Reference: PTXBC063822

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AGRINL; AGRNL; PIKA; pikachurin; EGF-like, fibronectin type-III and laminin G-like domain-containing protein; agrin-like protein; EGF like, fibronectin type III and laminin G domains		

For research use only.

GENE DESCRIPTION	EGF-like, fibronectin type III and laminin G domains
GENE ID	133584
ORIGIN SPECIES	Human
ACCESSIONS	BC063822
VECTOR	pDONR223
NCBI SYMBOL	EGFLAM
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr