

## C5orf28-chromosome 5 open reading frame 28 Gene

Reference: PTXBC063832

|                               |                                                                   |                     |                |
|-------------------------------|-------------------------------------------------------------------|---------------------|----------------|
| Product type:                 | DNA                                                               | Storage Condition:  | -20°C          |
| Size/Quantity:                | 2 ug                                                              | Delivery Condition: | +4°C, Blue Ice |
| Synonyms                      | transmembrane protein C5orf28; C5orf28; transmembrane protein 267 |                     |                |
| <b>For research use only.</b> |                                                                   |                     |                |

---

|                         |                                                                                                   |
|-------------------------|---------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | chromosome 5 open reading frame 28                                                                |
| <b>GENE ID</b>          | 64417                                                                                             |
| <b>ORIGIN SPECIES</b>   | Human                                                                                             |
| <b>ACCESSIONS</b>       | BC063832                                                                                          |
| <b>VECTOR</b>           | pDONR223                                                                                          |
| <b>NCBI SYMBOL</b>      | C5orf28                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                    |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon )

atggcatccgagactgaaaagaccatgctttactgcagactgttagcactgaatctcttattccagccttggctctggggcattttgcctctgtagctgacagacttctcagtttccacaattcagcaaatgactggctctgctctcagataatgcagtacattgttaattggcatgtggt  
catggcggttagtcactggaatcaagaagaagactgactttggagaatcatttttagctggatttttagcctctgttattgatgtagaccatttttctagctggatccatgctttaaggctgctttgactctccgagacatttcttcaactgttctactgtgattccggtgtgtgtctgacct  
gaaatttactatgcaccttttcaagctcaagactcatgggtcttctccctggatgttatttatactctggacttcacatcataccgagatgggattcgtcatggtttgtggatgccatttggaaaaactctcttggcattctggtttatgtaataatcacatcatctttacctcacatctgttcat  
tcgttatgtatttaacaggaaaccagacaaatgatgtcttcaaaacatggagttcgtattgatgtctga