

C1orf109-chromosome 1 open reading frame 109 Gene

Reference: PTXBC063843

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	uncharacterized protein C1orf109; chromosome 1 open reading frame 109		
For research use only.			

GENE DESCRIPTION	chromosome 1 open reading frame 109
GENE ID	54955
ORIGIN SPECIES	Human
ACCESSIONS	BC063843
VECTOR	pDONR223
NCBI SYMBOL	C1orf109
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgactcaagaccggcctctgttgcctgcaggaggcgtgaagaagtgttcccggtgtggaggagcagcaggccctgtggcagagtgccttgcgggactgccagccctctgtctccctcagcaacctggcggaacagctgcaggccgcacagaacctgcggtttgaggat
gtgccggcgcttcggcctccagatttaaaagagcggctgaggcgtgaagcagctgggtggctgtgacatcgtctggacaagctaggggaaaggctagccatcctcctcaagggtgcgagacatgttcagcagccatgtggagcaggtttcagatctatgcaacacgcagacac
agttggcattgatgtctgtcctgcagccttcagcagtgagccctctgtggctgacatgttggatgttgcaggatattgagagacattatcgaaagtcgtacctgaagagaaagtatcttcttctgtctatccagtgaggagacttggcaaacatacaagcttgcacaggcctgggaccgaa
tttcaaaagacgaaccaagatctgtacaagatatcctattgaatgttccctctcctggaagagtaa