

## RANBP9-RAN binding protein 9 Gene

Reference: PTXBC063849

|                |                                                                                                                                                                                                                                 |                     |                |
|----------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------------|----------------|
| Product type:  | DNA                                                                                                                                                                                                                             | Storage Condition:  | -20°C          |
| Size/Quantity: | 2 ug                                                                                                                                                                                                                            | Delivery Condition: | +4°C, Blue Ice |
| Synonyms       | BPM-L; BPM90; RANBPM; RanBP7; ran-binding protein 9; Ran Binding Protein in the Microtubule organizing center; novel centrosomal protein RanBPM; ran binding protein, centrosomal; ran-binding protein M; RAN binding protein 9 |                     |                |

**For research use only.**

|                         |                                                                                                  |
|-------------------------|--------------------------------------------------------------------------------------------------|
| <b>GENE DESCRIPTION</b> | RAN binding protein 9                                                                            |
| <b>GENE ID</b>          | 10048                                                                                            |
| <b>ORIGIN SPECIES</b>   | Human                                                                                            |
| <b>ACCESSIONS</b>       | BC063849                                                                                         |
| <b>VECTOR</b>           | pDONR223                                                                                         |
| <b>NCBI SYMBOL</b>      | RANBP9                                                                                           |
| <b>LENGTH</b>           | ProteoGenix bp                                                                                   |
| <b>SEQUENCE PRIMERS</b> | Forward primer M13R (5'→3'): GTAAACGACGGCCAGT<br>Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG |

### ORF (Only for reference, product provided without STOP codon)

atcgaggagtggaacacaaatccaggcacagatagatcgatttctatcgagatcgagaaggagaatggcagaccatgatacaaaaatggttcatttatttagtccaccatgggtactgtgccacagcagaggccttggcagatcacagaccagaccgttctagaagaattagct  
tccattaagaatagacaagaattcagaattggtattagcaggaagaatgggagaagccattgaaacaacacacagttatacccaagtttacttgaagaaatcctaatctctttcacattaaaagtcgctcagtttatagaatggtgaatggtacagatagtgagtagctgttgggag  
gccgaagtcacaaagtcacagacagttatctctgtagtccctgaccttttagttagtccagtagtagagcccccagccatggaatgaatatacccaatttagcatcaggcaaggaagcaccgcacattttcaggttttgaagttgtagtaattggtgaatcaataaagcacatcaatcatattgc  
catagtaataaacaccagtcacacacttgaatgtaccagaactaaacagtataaataatgtcaagatcacagcaagttataaacttaccagtaattgatgtagacatggaacagatcactcctcaatggagttggagaacttcatcaatgggttctcaatgggtgctctaaacatgaccac  
gaaatggaagattgtgacaccgaaatggaagtgtattcaagtcagttgagacgccagttgtgtggaggaagtcaggccgcatagaagaatgatccacttggagcagagctgcaagcaatgagtgaaacagtaaggagagactgtggcaagaacactgcaacaaaaaatgttga  
ggatgcattcagtcactagcatattcagatccctggaacagccagttggaatcagcttgaccgattcagagagaacctgtgtgctcagctcttaacagtgcattattagaacccacaatctgccaaagcaacctccacttgccttagcaatgggacagccacacaaatgtctaggact  
gatggctcagatcaggaattggatcctgcgcatgtgccagtggaagactacacattag