

KIAA0892-KIAA0892 Gene

Reference: PTXBC063863

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KIAA0892; MAU2L; SCC4; mau-2; MAU2 chromatid cohesion factor homolog; cohesin loading complex subunit SCC4 homolog; protein MAU-2; sister chromatid cohesion 4; MAU2 sister chromatid cohesion factor		

For research use only.

GENE DESCRIPTION	KIAA0892
GENE ID	23383
ORIGIN SPECIES	Human
ACCESSIONS	BC063863
VECTOR	pDONR223
NCBI SYMBOL	KIAA0892
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggacaacgcggaagcccagttcaccacggccctgcggctaccaaccaccaggagctgtgggcttcacgtcaccaacctggcgagtggtatatacgggaaggaatagacaccaagagctctacagctctgctggagaggatcaacccggaccacagcttccctgtcagctcgc
actgctccgagcagccgcttctatgtgctgtggcttctctcttccaggagcgtacaacgaggccaagcgatttctcgggaaactctgaagatgtccaatgctgaggacctgaacggctcacagctgctccctgtgcttctgggacacattctatgtgctggaaaccaca
gggagagtacaacatgggtgctgctgcatgcagctgcagcagaagatcccgacatgtcgttacagctgtggtcgtcagcactgctgagagacctgaataaagcctgtgggaacgcatggatgccatgaagccgcccagatgccagaacttctgcagcagctgctccagga
ccacattgaggcctgcagctccccgaacacaacctcatcagctggacagacggtccacccccgtgcagttccaagctcagaatggaccaacaccagcctggcagcctctgtga