

C6orf182-chromosome 6 open reading frame 182 Gene

Reference: PTXBC064365

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	C6orf182; bA487F23.2; cep57R; centrosomal protein CEP57L1; centrosomal protein 57kDa-like 1; centrosomal protein 57kDa-like protein 1; centrosomal protein of 57 kDa-related protein; cep57-related protein; centrosomal protein 57 like 1		

For research use only.

GENE DESCRIPTION	chromosome 6 open reading frame 182
GENE ID	285753
ORIGIN SPECIES	Human
ACCESSIONS	BC064365
VECTOR	pDONR223
NCBI SYMBOL	C6orf182
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggattctgaattaatgcatagtatagtaggaagctatcataaacctccagaagagtagttgttcctcattcaccagagaatgaacctctcagaattgccatctgcaacctagaagttacctctcctaaagatactcatagcccaaatagccaagctcttatttagccttaaaactcttcagg
aaaaaattcatcgttttagagctggagagaacacaagctgaagataaacctgaacattcttccagagaagcagcacagtagataagaaggccttagagaatgaacaaatgagagaaatctggcacaccaggagctgataaagcagaaaaagataagtatagcttaagctcagccagtc
tcgttgactcttctagagaagcaactagaatatacaagagaatggttctcaacgtagagcgagaaaaagaacatgatcctagaacaacaggccagcttcagagggaaaaagaacaagatcagatgaagctgtatgcaaaacttgaaaagcttgatgtcttagaaaaagagtggttcagact
tacaacaactcagaaaaactgctgaggacaagattaaacatttagaagaaaaacttaaggagaagaacatcagcgttaagctatttcaagacaagctctgagcttcaactggacttgaaatcagtaaaattataatgtcttcagtttcaaatctaaagcactccaaggaaaaagaagaatcttc
aaagggtgtgcataata