



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

WISP2-WNT1 inducible signaling pathway protein 2 Gene

Reference: PTXBC064379

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CCN5; CT58; CTGF-L; WNT1-inducible-signaling pathway protein 2; CCN family member 5; connective tissue growth factor-like protein; connective tissue growth factor-related protein 58; WNT1 inducible signaling pathway protein 2		

For research use only.

GENE DESCRIPTION	WNT1 inducible signaling pathway protein 2
GENE ID	8839
ORIGIN SPECIES	Human
ACCESSIONS	BC064379
VECTOR	pDONR223
NCBI SYMBOL	WISP2
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr