



Reference: PTXBC064383

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	zinc transporter SLC39A9; ZIP-9; ZIP9; zinc transporter ZIP9; solute carrier family 39 (zinc transporter), member 9; zrt- and Irt-like protein 9; solute carrier family 39 member 9		

For research use only.

GENE DESCRIPTION	solute carrier family 39 (zinc transporter), member 9
GENE ID	55334
ORIGIN SPECIES	Human
ACCESSIONS	BC064383
VECTOR	pDONR223
NCBI SYMBOL	SLC39A9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]