

SCLT1-sodium channel and clathrin linker 1 Gene

Reference: PTXBC064428

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CAP-1A; CAP1A; sodium channel and clathrin linker 1; sodium channel-associated protein 1		

For research use only.

GENE DESCRIPTION	sodium channel and clathrin linker 1
GENE ID	132320
ORIGIN SPECIES	Human
ACCESSIONS	BC064428
VECTOR	pDONR223
NCBI SYMBOL	SCLT1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgcagaaatcgactttctgagagagcaaatcgaagactaaatgaagatttaggcggtatcaaatgaaagttttccaaatattcatctgtacagaaagctgtctgccaaggagaaggagacacatttgaaaacctagatttgaccaaagcttttagctcctctgttactgagta
tgataaacacctaggagaactaaatggcgagctgaaatattaccaggcataa