



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

LY9-lymphocyte antigen 9 Gene

Reference: PTXBC064485

Product type:

DNA

Storage Condition:

-20°C

Size/Quantity:

2 ug

Delivery Condition:

+4°C, Blue Ice

Synonyms

CD229; SLAMF3; hly9; mLY9; T-lymphocyte surface antigen Ly-9; SLAM family member 3; cell surface molecule Ly-9; signaling lymphocytic activation molecule 3; lymphocyte antigen 9

For research use only.

GENE DESCRIPTION	lymphocyte antigen 9
GENE ID	4063
ORIGIN SPECIES	Human
ACCESSIONS	BC064485
VECTOR	pDONR223
NCBI SYMBOL	LY9
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggtggcaccacaagagtcacacagatgactgggctctggcccttctccagtaagccacagaggagtcagctgcgaataatctctctgtctctacagacctctctctctctctctatgggactaagagcctctggaaaggactcagcccccaacgtggtgtcagggtatcatgggggttccgtgactctccccataacatctcagttagacacagagattggaacgtctcatctggatgggtcccaaaaatgctctgtcttgcagctcccaaaagaaatgtaaccattatgtgcaaaagctacctggggccgactagacatcaccaggtgagattactccctgtcatcagaacatctgactctgaatgatgcaggatctctcaaaagccagataaaccaaaaggaatttgaagctcaccatgaggaggaaattccacctgttctgctatgcacacattattgaaaagtgttcgcagctctcaggggtgaccaccgcacactcctgggggacgcccgtgttcattcatcagcaccttgc

More info: Please contact us for any question: contact@proteogenix.fr

Page 1/1