

CORO6-coronin 6 Gene

Reference: PTXBC064514

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	coronin-6; clipin-E; coronin, actin binding protein 6; coronin-like protein E; coronin 6		
For research use only.			

GENE DESCRIPTION	coronin 6
GENE ID	84940
ORIGIN SPECIES	Human
ACCESSIONS	BC064514
VECTOR	pDONR223
NCBI SYMBOL	CORO6
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgccgagtccttgggggtgttctgtgactgcatgcccgcgcagcgggaacaacttcgaggagccagtggcactgcaggagatggacacaagcaacggggctctattgccctttacgatcccgactccagcatcgtctacgtgtggcaagggcgacagcagcattcgtacttt
gagattaccgacgagccgcttctgtgactacgtgaacacgttcagcagcaaaagagccgcagcggggcatgggtttcatgccccaaaggggactggatgtcagcaagtgtgagatcggcgttctacaagctacacgaagaaagtgtgaacctatcatgactgtccccgaa
gtcagacctcttccaggacgatctgtaccggatagccagggccggagccggccctagaagcggagcaatggctatccggccaggacgccgaacccgtgtctattcgtgaggacggctatgtcccccaagcaccgcgagctccgggtcacgaagcgcaacatcctggacg
tgcggccgctccggccccgcccgcagccagtcggccagcagcccccttgcgagcacaccctggagacgctgctggaagagatcaaggccctccgagcgggtgcaggccagcagcgcacacggctctggagaacatgctgtgcgagctggtggacggcagc
gactag