

TMEM53-transmembrane protein 53 Gene

Reference: PTXBC064520

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NET4; transmembrane protein 53; novel DUF829 domain-containing protein		
For research use only.			

GENE DESCRIPTION	transmembrane protein 53
GENE ID	79639
ORIGIN SPECIES	Human
ACCESSIONS	BC064520
VECTOR	pDONR223
NCBI SYMBOL	TMEM53
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcctggcagagctggactacacatcgagatccggatcagccctgctggagccagaagaacagccccagccaggtgggaaggaggcagaactcggcagcctgtggtgattctcttggcctggggtggcctgcaaggacaagaacctggcaagtacagtccatctaccac
aaaaggggctgcatcgtaatccgatacacagccccgtggcacatggctcttctccgagtcactgggtatcccttcacttcgtgttttggccagaagctgctcagctcttctgattatgagattgagaaggagccccctgctcttccatgtcttcagcaacgggtgctatgctgtaccgt
acgtgctggagctcctgagacccgtcgttctgcccctgctgtgtggggacacatcttgcacagcgtcctctgtgacagcaacctggttagggcctctcggggccctggcagccatcctggagcggccggcccatgctgcgcctgttgcctgtgtggccttgcctgtgtgtcgt
cctgttcacgtcctgttgcctccatcacagcccttccacacccacttctatgacaggclacaggacgggctcgtcgtggcccgagctctactctactcggaggctgacgaagtagtctcggccagagacatagaacgcatgtgtggaggcacgcctggcagccgggtcctggc
gcgttctgtggttctgtatctgcacacgtcagccacctcgtgactaccctacttactacacaagcctctgtgtcgaactcgtccgtcgtga