

DYNC1H1-dynein, cytoplasmic 1, heavy chain 1 Gene

Reference: PTXBC064521

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CMT2O; DHC1; DHC1a; DNCH1; DNCL; DNECL; DYHC; Dnchc1; HL-3; SMALED1; p22; cytoplasmic dynein 1 heavy chain 1; dynein heavy chain, cytosolic; dynein, cytoplasmic, heavy polypeptide 1; dynein cytoplasmic 1 heavy chain 1		

For research use only.

GENE DESCRIPTION	dynein, cytoplasmic 1, heavy chain 1
GENE ID	1778
ORIGIN SPECIES	Human
ACCESSIONS	BC064521
VECTOR	pDONR223
NCBI SYMBOL	DYNC1H1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgatcagtaaaatgctgaagatgcagatgttgaggatgaggacacctggcctacgcagagaccgagaagaagacgaggacagactccacgtccgacggcgccctgctggatgctggacactgcacaccaccgctccaactggctgacactatccccagacgctgagcca
cctcaagcgaccgtggagaatatcaaggatccttfttcaggtttttgagagagaagtgaagatggcgcaaaagctgtctcaggacgttcgccaggacctgcagatgtctgctcaggtgtgcgaaggaaaagaagaagcagaccaactacttgcgcacgtgatcaacgactagttaa
agggatcttgcctcggagctgttccactacacggctgctcggcatgaccgtatccagttgggtgtccgacttcagcagagagatcaaacagctgcagaacatctcactggcagctgcatctggtggcgcaaggagctaaaggtgaaggcctcctgacgagctctgggtggtcag
cagctgctcgggtgggtggagtgctctgggaaaaaacagggccaggtctga