

FLJ35220-hypothetical protein FLJ35220 Gene

Reference: PTXBC064545

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	endonuclease V; hEndoV; inosine-specific endoribonuclease		

For research use only.

GENE DESCRIPTION	hypothetical protein FLJ35220
GENE ID	284131
ORIGIN SPECIES	Human
ACCESSIONS	BC064545
VECTOR	pDONR223
NCBI SYMBOL	FLJ35220
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggccctggagcgggcgaggccgagggaacgctgtcactgtggaacgggagcaagctcggctgaaggccacgtcgtagaccgggacaccgaggcgtggcagcgagaccccgcttctcgggtcagaggggtcggggcggtgacgtgtcctcgtgaaagg
ggacagtggtccgcttgtgtccctgtgtgtcagcttccctgagctcaggtctcttctgtggtggaacggggtactccaccaccgaggcttgggggtggcctgccaccttggcgtcttacagacctgccgtgtgttgggggtggccaaagaaactctcaggtggatgggtgga
gaacaacgctgcacaaggagaagatccgactcctcagactcagggagactcattccctctgtctgggagactctgggactgtcctgggaatggcctgaggagccacgaccgcagcaccagggccctctacatctcgttggccacaggatgagcctggaggccgtgtgcgcct
gacttctgtctcaggttccggatcccagagccggtgcgcccaggtgacatctctccgagagcacatccgcaagtcgtgggactccccgggcccacccacccgaggagcccgaaggcgagagggcagtggtatgccccaaaggagactccggagagtcctcaggtgag
ggccagccccacagaccacagcccagggccagggcagcccaaggccaggtccagagcagggcggaaggactggcagtag