



Antibodies and Proteins

ProteoGenix SAS
Espace Européen de l'Entreprise
15 rue de La Haye
67300 Schiltigheim, France
(0)3 90 20 54 70 - Fax +33 (0)3 88
contact@proteogenix.fr
www.ProteoGenix.science

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88
contact@proteogenix.fr
www.ProteoGenix.science

PRH1-proline-rich protein HaeIII subfamily 1 Gene

Reference: PTXBC064553

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	Db-s; PIF-S; PRH2; Pr1/Pr2; salivary acidic proline-rich phosphoprotein 1/2; PRP-1/PRP-2; parotid acidic protein; parotid double-band protein; parotid isoelectric focusing variant protein; parotid proline-rich protein 1/2; protein C; proline rich protein		

For research use only.

GENE DESCRIPTION	proline-rich protein HaeIII subfamily 1
GENE ID	5554
ORIGIN SPECIES	Human
ACCESSIONS	BC064553
VECTOR	pDONR223
NCBI SYMBOL	PRH1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcttctgattctgctgagctggccctctggccttcagctcagctcaggtatfaaataagaagatgtagccaggagaaatggttcctctgtaatacagatggagagactctgagcaggttctatagataggagcgtcaggagaccactcttggagagacgaactcaacccctctgctgtgtg
atgggaacacaggatgatggcctcagcagggaccaccccaacaaggaggccagcagcaacaaggttcaccactctctcagggaagccacaaggagagccagcagcaacaaggtccaccactctctcagggaagccacaaggaccaccccaacaggg
agagccatccctctctctcaaggagccacaaaggaccaccccaacaggagagccalccccgtctctctcagagaaagccacaaggaccaccccaacaggagagccalcgaaaggtctctccccactctctcagggaagcccaaggaccactccccaaaggggcccgcc
acaagagctctccacagggggcagctctctcagtaa

More info: Please contact us for any question: contact@proteogenix.fr