

ELL-elongation factor RNA polymerase II Gene

Reference: PTXBC064558

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	ELL gene (11-19 lysine-rich leukemia gene); RNA polymerase II elongation factor ELL; C19orf17; ELL1; MEN; PPP1R68; eleven-nineteen lysine-rich leukemia protein; elongation factor RNA polymerase II; protein phosphatase 1, regulatory subunit 68; elongation		

For research use only.

GENE DESCRIPTION	elongation factor RNA polymerase II
GENE ID	8178
ORIGIN SPECIES	Human
ACCESSIONS	BC064558
VECTOR	pDONR223
NCBI SYMBOL	ELL
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgacttggaaggtgtgtgggggtggatggcgctgtgactcaggcccgatcacatgggggtcactgctccatcagcagtcaggcaggagacagaaatgcaaaatacagatgtctattttctaaactgggggtggaggggtgcctcctgacagcttcaagagctga
ctgtggagggccacagtgatccaccctgcccagagctctctcttgggaccccggttccagagctccgaggaaactga