



Reference: PTXBC064567

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	GE1; Ge-1; HEDL5; HEDLS; RCD-8; RCD8; enhancer of mRNA-decapping protein 4; autoantigen Ge-1; autoantigen RCD-8; human enhancer of decapping large subunit; enhancer of mRNA decapping 4		

For research use only.

GENE DESCRIPTION	enhancer of mRNA decapping 4
GENE ID	23644
ORIGIN SPECIES	Human
ACCESSIONS	BC064567
VECTOR	pDONR223
NCBI SYMBOL	EDC4
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]