



Reference: PTXBC064581

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	RING finger protein 214		
For research use only.			

GENE DESCRIPTION	ring finger protein 214
GENE ID	257160
ORIGIN SPECIES	Human
ACCESSIONS	BC064581
VECTOR	pDONR223
NCBI SYMBOL	RNF214
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]