

KIAA1984-KIAA1984 Gene

Reference: PTXBC064582

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KIAA1984; PARF; bA216L13.7; coiled-coil domain-containing protein 183; coiled-coil domain containing 183		

For research use only.

GENE DESCRIPTION	KIAA1984
GENE ID	84960
ORIGIN SPECIES	Human
ACCESSIONS	BC064582
VECTOR	pDONR223
NCBI SYMBOL	KIAA1984
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggtggtgcgggagaagctgcgaagtacgtcttcgaccgctgaacatgcacaacctactgatccacctggtgcggcgcgcggggcagaagctggagagcatgcagctggagctggacagctgcggagccagcccgacgccagcaaggaggagctgcggctgctgcag
atcatccgccagctggagaacaacatcgagaagacaatgatcaagatcatcaccagccagaacatccacctgctgatttggacctgctggattatctgaagacagtgctggcaggatccccattgagctggacaagctgcagaacctctgtgtcaactactgtcagagctgtcgatat
gaagatcatgtcccaagatgccatgatgatcacggatgaggtcaaggaaacatgaggcaaggaggcgtccttcacgaggagcgcggggcaaggagaaacggctcaaccagcagaagaagctgatcgacaagatccacacgaaggagaccagcgagaagtaccgccggggc
cagatggacttggacttcccctgaaacctgatgagcacggagacctgaaattgaggagaaaagagacctccacagcagaatggaataaccagtcgggcgtgactgctgtgtgtggagaaggtcaagagtgctgtacggctgtcactctctggacatcactag