

YIPF1-Yip1 domain family, member 1 Gene

Reference: PTXBC064590

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein YIPF1; DJ167A19.1; FinGER1; YIP1 family member 1; Yip1 domain family member 1		

For research use only.

GENE DESCRIPTION	Yip1 domain family, member 1
GENE ID	54432
ORIGIN SPECIES	Human
ACCESSIONS	BC064590
VECTOR	pDONR223
NCBI SYMBOL	YIPF1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcagcagtagatgacttgcatttgaagaatttggcaatgcagccacttctctgacagcaaacccagatgccaccacagtaaacattgaggatcctgtgaaaccccaaacatcagccaggatcccaagaggctcaggaagagaagaagatgatgagttactgggaatgatgac
tctgacaaaactgagttacttctgctgacagagaagaaagctccccctctggacatttgaatactaccaaacattcttggatgtgacacctaccaggctcttggacagaattaaaggatctcttggcaataccgggaaaaacttggagggttatataccgcagcaatccagatctctatggccc
ctttggatatgtgccaggttggctcttggccatagcaattagtggaatcttccaacttcttgatccatctgggagagaagacgtaccattatgtgcccgaattccgaaaaagtgtccatagcagctaccatcatctatgccatgacctggctgttccctctgacactctggggttccatgtggaga
aacagcaaaagtatgaacatgctctctatttcttggagattgtgtgtctatggatctccctcttcatctatccccaccgcaatactgtgattatccccagaaagctgttgcgttgattctatgcatgattgcctgggcatctcaggatctcttggcaatgacatttggccagctgttc
gtgaggataaccgacgcttgcattggccacaattgtgacaattgtgtgctccatgtctgttctgtgggctgttgcatactttttgatgcaccagagatggaccatctcccaacaactacagctactccaaaccaaacagttgtgcagccaagtcagctaa