

VPS37D-vacuolar protein sorting 37 homolog D (*S. cerevisiae*) Gene

Reference: PTXBC064621

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	VPS37D, ESCRT-I subunit; ESCRT-I complex subunit VPS37D; WBSCR24; vacuolar protein sorting-associated protein 37D; Williams Beuren syndrome chromosome region 24; Williams-Beuren syndrome critical region protein 24; vacuolar protein sorting 37 homolog D; v		

For research use only.

GENE DESCRIPTION	vacuolar protein sorting 37 homolog D (<i>S. cerevisiae</i>)
GENE ID	155382
ORIGIN SPECIES	Human
ACCESSIONS	BC064621
VECTOR	pDONR223
NCBI SYMBOL	VPS37D
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgcatcgtggagtcctccactgcgcgtgggtgctgctgaggtgagctagaagaggcggagcaggagcagagagcagatggagcagctgctgctggggagcaaaagcctggaggccttctgctgcttccagcgtggccgcgcctggccacctgaggcggacgcag
gcagagaagctgcaggagctgctgctggcgtcgggagcgttctgcccagccggccccacctcggctgctgatcccccaaatcttcccggtgcagctgtcttgcctggggccgcccggggccaccagcagtgccccggagcctgcccccttggactccgcccagtgcc
ccactgaagggtccccgggtgccccctcggccggccccctgctgagccctcggccctcgagccagagccccccaccggtag