



Reference: PTXBC064663

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	serine/threonine-protein kinase NLK; nemo like kinase		
For research use only.			

GENE DESCRIPTION	nemo-like kinase
GENE ID	51701
ORIGIN SPECIES	Human
ACCESSIONS	BC064663
VECTOR	pDONR223
NCBI SYMBOL	NLK
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]