

TNFAIP3-tumor necrosis factor, alpha-induced protein 3 Gene

Reference: PTXBC064689

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	AISBL; OTUD7C; TNFAIP2; tumor necrosis factor alpha-induced protein 3; OTU domain-containing protein 7C; tumor necrosis factor inducible protein A20; tumor necrosis factor, alpha induced protein 3; zinc finger protein A20; TNF alpha induced protein 3		

For research use only.

GENE DESCRIPTION	tumor necrosis factor, alpha-induced protein 3
GENE ID	7128
ORIGIN SPECIES	Human
ACCESSIONS	BC064689
VECTOR	pDONR223
NCBI SYMBOL	TNFAIP3
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggctgaacaagtcctcaggctttgtattgagcaatgacggaaagctgtgaagatcgggagagaaactcagaagacattttaacctactaatgggacattcatcttttaaacatgcaccgatacacactggaatgtcagaactgccagtttgcctcagtttcgggagat
catccacaagccctcatcgacagaacatccaggccaccctgaaagccagaagaactcaactggtgtcgagaagtcggaaagcttggcgctgaaacgaacgggtgacggcaattgcctcatgcatgccacttcagtcacatgtggggcgttcaggacacagactggtactga
ggaaggcgtgttcagcagctcaaggaaacagacacagcaactttaattccgctggcaactggagctctcaaatctcaggaattgttgaacggggccttgcctatgatactcggaactggaatgatgaatgggacaatcttcaaaatggcttccacagacacccatggccgaa
gtggacttcagtacaactcactggaagaatacacatattgtccttgcacatcctcagaagccaatcattgtcattcaggtgagatgcctgcagatcacggatctgtacttaaatgcttcagccttatgccttggctcctggagaaaaccacactgccaaagtcaggtacagagttcaat
ggaatttga