



Reference: PTXBC064694

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	hWNT5A; protein Wnt-5a; WNT-5A protein; wingless-type MMTV integration site family, member 5A; Wnt family member 5A		

For research use only.

GENE DESCRIPTION	wingless-type MMTV integration site family, member 5A
GENE ID	7474
ORIGIN SPECIES	Human
ACCESSIONS	BC064694
VECTOR	pDONR223
NCBI SYMBOL	WNT5A
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]

More info: Please contact us for any question: contact@proteogenix.fr