



Antibodies and Proteins

ProteoGenix SAS

Espace Européen de l'Entreprise

15 rue de La Haye

67300 Schiltigheim, France

Tel. +33 (0)3 90 20 54 70 - Fax +33 (0)3 88 56 16 88

contact@proteogenix.fr

www.ProteoGenix.science

KIAA0247-KIAA0247 Gene

Reference: PTXBC064697

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	KIAA0247; DRAGO; sushi domain-containing protein 6; drug-activated gene overexpressed protein; sushi domain containing 6		

For research use only.

GENE DESCRIPTION	KIAA0247
GENE ID	9766
ORIGIN SPECIES	Human
ACCESSIONS	BC064697
VECTOR	pDONR223
NCBI SYMBOL	KIAA0247
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

[illegible]