



Reference: PTXBC064803

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	MMP-14; MMP-X1; MT-MMP; MT-MMP 1; MT1MMP; MTMMP1; WNCHRS; matrix metalloproteinase-14; matrix metalloproteinase 14 (membrane-inserted); membrane type 1 metalloproteinase; membrane-type-1 matrix metalloproteinase; matrix metalloproteinase 14		

GENE DESCRIPTION	matrix metallopeptidase 14 (membrane-inserted)
GENE ID	4323
ORIGIN SPECIES	Human
ACCESSIONS	BC064803
VECTOR	pDONR223
NCBI SYMBOL	MMP14
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]