



Reference: PTXBC064806

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	CD335; LY94; NK-p46; NKP46; natural cytotoxicity triggering receptor 1; NK cell-activating receptor; lymphocyte antigen 94 homolog (activating NK-receptor; NK-p46); natural killer cell p46-related protein		

For research use only.

GENE DESCRIPTION	natural cytotoxicity triggering receptor 1
GENE ID	9437
ORIGIN SPECIES	Human
ACCESSIONS	BC064806
VECTOR	pDONR223
NCBI SYMBOL	NCR1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]