



Reference: PTXBC064807

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	PTPSTEP; STEP; STEP61; tyrosine-protein phosphatase non-receptor type 5; neural-specific protein-tyrosine phosphatase; protein tyrosine phosphatase, non-receptor type 5 (striatum-enriched); protein-tyrosine phosphatase striatum-enriched; striatal-enriched		

For research use only.

GENE DESCRIPTION	protein tyrosine phosphatase, non-receptor type 5 (striatum-enriched)
GENE ID	84867
ORIGIN SPECIES	Human
ACCESSIONS	BC064807
VECTOR	pDONR223
NCBI SYMBOL	PTPN5
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

[illegible]