

C3orf33-chromosome 3 open reading frame 33 Gene

Reference: PTXBC064898

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	protein C3orf33; AC3-33; AP-1 activity suppressor; chromosome 3 open reading frame 33		

For research use only.

GENE DESCRIPTION	chromosome 3 open reading frame 33
GENE ID	285315
ORIGIN SPECIES	Human
ACCESSIONS	BC064898
VECTOR	pDONR223
NCBI SYMBOL	C3orf33
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggcggggcagcccgccgaccggctcgcgtctgcccagacaaggacgggaatggagcccaactcgtggctcggatctcgcagtgggcagacgaccactcgcgctagtcgggaacatcagcactggaatggccatagctggaataatgttacttttgaagcattcgactgacat
caaaattacaagctcttcggatattccagtagaatttataagaagaatgttaactacgtggagattacgcccgaataactgagaatgggttagaattgaacatatactatttactattatagcttcattgagaaaagagccacgtggtgctttgctggttaagtggctggagtagaact
cgctgaaactgggaagcattggttcaaaaaagagctaaaacctcccaacttactgttccaaactcttggaaaggagaattcagcactctttgctatctctgtggaataagggtggatatttcagcgtgaatctgaatgaagaaatttgagaagagggccttggcaaacgttctgtttaaag
ggcttaaatgatctaaatctactggacagttcacagaacttacttaagctgaattaacagccttaaaaaaggagaagggaatgatgaaggagagactctgaaaaagaagttacttagaaaaattcaagattcctggagagaatatgaaaaaggacagtttttttaaaacaacagg
atcagatttcagcttgaaaaaagaaagttattatgaaaacttaaaaggacttatgaaatatggaagacaacatgaacaactgctccttaactgaagttcagagaactataagtcgataaaacttcgtagaaaaagggtga