

SCYL1BP1-SCY1-like 1 binding protein 1 Gene

Reference: PTXBC064945

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	SCYL1BP1; NTKLBP1; RAB6-interacting golgin; N-terminal kinase-like-binding protein 1; NTKL-binding protein 1; SCY1-like 1-binding protein 1; SCYL1-BP1; SCYL1-binding protein 1; hNTKL-BP1; golgin, RAB6 interacting		

For research use only.

GENE DESCRIPTION	SCY1-like 1 binding protein 1
GENE ID	92344
ORIGIN SPECIES	Human
ACCESSIONS	BC064945
VECTOR	pDONR223
NCBI SYMBOL	SCYL1BP1
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

```
atgagctgggcagcagtggtggcagtcgcggctcgcgagattgggcactttgggggtgccggcgccggccgatggcgcaaggtgggcaggtctctgagagggaactgaggagactaaagcagactaaagatccattgaaccacagcagctctccccgcgaagaaaagtc
gacaacaactcagcagaaaaagccctgtgagagcaaaagccaaaaactgggctcaagatggatcaactcattctccagagcagctcttcagcaccacaaacagagagtttaacgttcaaaaaccacctttctccctactctccgagtcattcactctcactccccgttggtg
atggacaaccacagggcattgaaagtcagccaaaggaactgggacttgagaattccatgatggtcacaacaatgttgagattctacctcaaaagccagattgcaaatggagaaaaaagaaatgggaattgcaagaaaaatctcgttgggaagtcctccaacaagaacggctaagg
agagaaaaataaacgtaaaaaagctctttggcctaagctattgcgaaagatctaaagaactcagccagagaccatgaaactaaagcggatccagaaggagttgcaggcttagatgacatggtgtcagctgacattggaattctcaggaaccggatgatcaggccagcttagactatc
atagctcgggga
```