

TMEM205-transmembrane protein 205 Gene

Reference: PTXBC064948

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	UNQ501; transmembrane protein 205; MBC3205		
For research use only.			

GENE DESCRIPTION	transmembrane protein 205
GENE ID	374882
ORIGIN SPECIES	Human
ACCESSIONS	BC064948
VECTOR	pDONR223
NCBI SYMBOL	TMEM205
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atggaggaaggcgggaacctaggagcctgattaagatgggtccatctactgtgtcaggtgcctggggcatgcaaatgtgggtgaccttgcctcaggcttctgcttttccgaagccttccccacatacttggactagtcagagcaaaactttcccttctactccatctccatg
ggctgtgcttcatcaacctctgcatcttggcttcacagcatgcttgggtcagctcacattctgggagggccagccagctttacgtgtgttctgagccttacgtggccactgtcaacgcccgtgctggaaacccgcaccacagctgcatgtggccctgcaaacgtggagaagga
gcgagggcctgggtggggaggtaccaggcagccaccagggtcccgaacctaccgccagctgcgagagaaggaccccaagtacagtgctctccgcagaatttctccgtaccatgggctgtctctcttggcaatctgggctgcgtcctgagcaatgggctgtctctgctggccttgc
ctggaaataaggagcctctag