

THOC7-THO complex 7 homolog (Drosophila) Gene

Reference: PTXBC065012

Product type:	DNA	Storage Condition:	-20°C
Size/Quantity:	2 ug	Delivery Condition:	+4°C, Blue Ice
Synonyms	NIF3L1BP1; fSAP24; hTREX30; THO complex subunit 7 homolog; NIF3L1-binding protein 1; Ngg1 interacting factor 3 like 1 binding protein 1; THO complex 7 homolog; functional spliceosome-associated protein 24; ngg1-interacting factor 3-like protein 1-binding		

For research use only.

GENE DESCRIPTION	THO complex 7 homolog (Drosophila)
GENE ID	80145
ORIGIN SPECIES	Human
ACCESSIONS	BC065012
VECTOR	pDONR223
NCBI SYMBOL	THOC7
LENGTH	ProteoGenix bp
SEQUENCE PRIMERS	Forward primer M13R (5'→3'): GTAAAACGACGGCCAGT Reverse primer T7P (5'→3'): TAATACGACTCACTATAGG

ORF (Only for reference, product provided without STOP codon)

atgggagccgtgactgacgacgaagttatacggagcgtctcctcatgatggagatgggtcctggagatgacggagaattaatctgctagtgaagagtttcattaaatggtgcaactctgggtcccagggaagagggatagccagtagcacaacgtatgctgagcagcgtgtctcaatgtgaa
ttttcaatgggcaaaactttactagtatgatgatgaatcagagaaatggaataatgaaaaaattacaaggaaatagaatgtagcatagctggagcagatgaaaaaattgctgagtgcaaaaagcaaatcttcaagcaaaacgaatagcaaaaaatcgccaagaatgatgctttggca
aaagtattcagcaccatccagacagggcatgagacattaaaggaaactagagcctctgggaaaaagaattagagcatcttcacacattaaagaagtggtgaagataagctggaattgagacggaacagtttcatgttctttagtaccatccatgaattcagcaaacattggaaaatgatga
aaaactctcagaggtagaagaagctcaggaagcaagcatggaacagatcctaagccatag